

Serum Plasma (heparin) 1~30°C

Vcheck Canine CRP 2.0

Dilute sample Use a 5 µl pipette to draw 5 µl of the sample and add into an assay diluent bottle.

Mix Close the bottle cap and shake for 5-6 times to mix thoroughly.

Measure Add 100 µl of the mixed sample to the sample well of the test device and press [START].

| | | |
|---------------|------------------|-----------------|
| Normal | Equivocal | Abnormal |
| < 20 mg/L | 20-30 mg/L | > 30 mg/L |

Serum Plasma (heparin) 1~30°C

Vcheck Feline SAA 3.0

Dilute sample Use a 5 µl pipette to draw 5 µl of the sample and add into an assay diluent tube.

Mix Close the tube cap and shake for 5-6 times to mix thoroughly.

Measure Add 100 µl of the mixed sample to the sample well of the test device and press [START].

| | | |
|---------------|------------------|-----------------|
| Normal | Equivocal | Abnormal |
| < 5 µg/ml | 5-10 µg/ml | > 10 µg/ml |

Serum 1~30°C

Vcheck cPL 2.0

Dilute sample Use a 25 µl pipette to draw 25 µl of the serum and add to the assay diluent tube.

Mix Use a 100 µl pipette to mix the sample with diluent by pipetting 5-6 times.

Measure Add 100 µl of the mixed sample to the sample well of the test device and press [START].

| | | |
|---------------|------------------|---------------------|
| Normal | Suspected | Pancreatitis |
| < 200 ng/ml | 200-400 ng/ml | > 400 ng/ml |

Serum Plasma (EDTA) 1~30°C

Vcheck fPL 2.0

Dilute sample Use a 25 µl pipette to draw 25 µl of the sample and add into an assay diluent tube.

Mix Use a 100 µl pipette to mix the sample with diluent by pipetting 5-6 times.

Measure Add 100 µl of the mixed sample to the sample well of the test device and press [START].

| | | |
|---------------|------------------|---------------------|
| Normal | Suspected | Pancreatitis |
| ≤ 3.5 ng/ml | 3.6-5.3 ng/ml | ≥ 5.4 ng/ml |

Serum (*DO NOT USE SST TUBE) 2~8°C

Vcheck T4

Dilute sample Use a 50 µl pipette to draw 50 µl of the serum and add to the assay diluent tube.

Mix & Incubate Use a tablet pipette to mix over 8 times until the tablet is **completely dissolved**. Incubate the sample mixture at room temperature for 10 minutes.

Measure Add 100 µl of incubated sample to the sample well of the test device and press [START].

| | | | | | |
|------------|-------------|-------------------|------------------|------------------------|--------------------|
| Dog | Low | Low normal | Normal | High | Therapeutic |
| | < 1.0 µg/dL | 1.0-2.0 | 1.0-4.0 | > 4 | 2.1-5.4 |
| Cat | Low | Normal | Gray zone | Hyperthyroidism | |
| | < 0.8 µg/dL | 0.8-4.7 | 2.3-4.7 | > 4.7 | |

Serum (*DO NOT USE SST TUBE) 2~8°C

Vcheck cTSH

Dilute sample Use a 100 µl pipette to draw 100 µl of the serum and add to the assay diluent tube.

Mix Use the same pipette to mix the sample with diluent by pipetting 5-6 times.

Measure Add 100 µl of the mixed sample to the sample well of the test device and press [START].

| | |
|---------------|-------------|
| Normal | High |
| < 0.5 ng/mL | ≥ 0.5 ng/mL |

Plasma (sodium citrate) 2~8°C

Vcheck D-dimer

Sodium citrate Use a 50 µl pipette to take 50 µl of sodium citrate and add it into a 1.5 ml tube.

Whole blood Add 450 µl of whole blood to the line of 0.5 ml. *Sodium citrate(50 µl):whole blood(450 µl)=1:9 ratio

Mix Mix the tube gently using a wrist snap making an 8-character shape.

Centrifuge Centrifuge at 3,000 rpm for 15 min.

Dilute Sample Use a 5 µl pipette to draw 5 µl of plasma and add it to the assay diluent tube.

| | |
|---------------|------------------------------------|
| Normal | Abnormal (TED/DIC probable) |
| < 0.3 µg/ml | ≥ 0.3 µg/ml |

Serum Plasma (heparin) 2~8°C

Vcheck SDMA

* Please place each sample and packaged kit component at room temperature (15~30°C) 30 minutes before use. Warming them up by hand or any other improper elevation of the temperature should be avoided as it may generate invalid results.

Sample Use a 100 µl pipette to draw 100 µl of the sample and add to a 1.5 ml tube.

Pretreatment Use a 25 µl pipette to draw 25 µl of the Pretreatment Buffer and add it to the 1.5 ml tube. * The end of the pipette tip should be applied onto the inner wall of the tube to prevent the pipette tip from clogging.

Mix Mix well using a vortex mixer. If not, tap with fingers 6-8 times. * Ensure that the sample and pretreatment buffer are well mixed to form white sediment.

Centrifuge Centrifuge at 6,000 x g for 5 min within 10 seconds from the mixing stage. * (10,000 RPM for 5 min with the centrifuge provided by Bionote)

Collect supernatant Collect 50 µl of the separated supernatant and transfer to a new 1.5 ml tube. * Be careful not to bring the sediments at the bottom.

| | | |
|---------------|-------------------|--------------------------------|
| Normal | Suspected | Kidney disease probable |
| ≤ 14 µg/dL | 14.1 - 19.9 µg/dL | ≥ 20 µg/dL |

Serum (*DO NOT USE SST TUBE) 2~8°C

Vcheck cCortisol

Dilute sample Use a 50 µl pipette to draw 50 µl of the serum and add to the assay diluent tube.

Mix & Incubate Use the same pipette to mix the sample with diluent by pipetting 5-6 times. Incubate the sample mixture at room temperature for 10 minutes.

Measure Add 100 µl of incubated sample to the sample well of the test device and press [START].

* Please refer to the product Instructions for Use to determine the results of ACTH ST, LDDST, or HDDST.

Serum (*DO NOT USE SST TUBE) 2~8°C

Vcheck cProgesterone

Dilute sample Use a 50 µl pipette to draw 50 µl of the serum and add to the assay diluent tube.




Mix Use the same pipette to mix the sample with diluent by pipetting 5-6 times.

Measure Add 100 µl of the mixed sample to the sample well of the test device and press [START].

| | | | | | |
|------------------------------|---------------------|-----------------|----------------------|------------------|-----------------------|
| Anestrus or proestrus | Pre-LH surge | LH surge | Post-LH surge | Ovulation | Post-ovulation |
| < 1.0 ng/mL | 1.0 - 1.99 | 2.0 - 2.99 | 3.0 - 4.99 | 5.0 - 12.0 | > 12.0 |

Serum 1~30°C




Vcheck Feline NT-proBNP

| Dilute sample | Mix | Measure |
|---|---|---|
|  |  |  |
| Use a 100 µl pipette to draw 100 µl of the serum and add to the assay diluent tube. | Use the same pipette to mix the sample with diluent by pipetting 5-6 times. | Add 100 µl of the mixed sample to the sample well of the test device and press [START]. |

| | |
|---------------|-----------------|
| Normal | Abnormal |
| < 100 pmol/L | ≥ 100 pmol/L |

Serum 2~8°C




Vcheck Canine NT-proBNP

| Dilute sample | Mix | Measure |
|---|---|---|
|  |  |  |
| Use a 100 µl pipette to draw 100 µl of the serum and add to the assay diluent tube. | Use the same pipette to mix the sample with diluent by pipetting 5-6 times. | Add 100 µl of the mixed sample to the sample well of the test device and press [START]. |

| | | |
|---------------|--------------------|-----------------|
| Normal | Suspected | Abnormal |
| < 900 pmol/L | 900 - 1,800 pmol/L | > 1,800 pmol/L |

Serum Plasma (heparin) 2~8°C




Vcheck Equine SAA

| Dilute sample | Mix | Measure |
|---|---|---|
|  |  |  |
| Use a 5 µl pipette to draw 5 µl of the sample and add into an assay diluent bottle. | Close the bottle cap and shake for 5-6 times to mix thoroughly. | Add 100 µl of the mixed sample to the sample well of the test device and press [START]. |

| | | |
|---------------|------------------|-----------------|
| Normal | Equivocal | Abnormal |
| < 10 mg/L | 10 - 20 mg/L | > 20 mg/L |

Serum 1~30°C




Vcheck Feline Tnl

| Dilute sample | Mix | Measure |
|---|--|---|
|  |  |  |
| Use a 100 µl pipette to draw 100 µl of the serum and add to the assay diluent tube. | Use the same pipette to mix the sample with diluent by pipetting 5-6 times. | Add 100 µl of the mixed sample to the sample well of the test device and press [START]. |

| | | |
|---------------|-------------------|-----------------|
| Normal | Suspected | Abnormal |
| < 0.18 ng/ml | 0.18 - 0.28 ng/ml | > 0.28 ng/ml |

Serum 1~30°C




Vcheck Canine Tnl

| Dilute sample | Mix | Measure |
|---|--|---|
|  |  |  |
| Use a 100 µl pipette to draw 100 µl of the serum and add to the assay diluent tube. | Use the same pipette to mix the sample with diluent by pipetting 5-6 times. | Add 100 µl of the mixed sample to the sample well of the test device and press [START]. |

| | | |
|---------------|------------------|-----------------|
| Normal | Suspected | Abnormal |
| < 0.1 ng/ml | 0.1 - 0.2 ng/ml | > 0.2 ng/ml |

Serum Plasma (heparin) 2~8°C






Vcheck eProgesterone

| Dilute sample | Mix | Measure |
|--|--|---|
|  |  |  |
| Use a 50 µl pipette to draw 50 µl of the sample and add to the assay diluent tube. | Use the same pipette to mix the sample with diluent by pipetting 5-6 times. | Add 100 µl of the mixed sample to the sample well of the test device and press [START]. |

| | |
|------------|-------------------------|
| Low | High (gestation) |
| ≤ 2 ng/ml | > 2 ng/ml |

2~30°C

Vcheck Ag, Ab

| Dilute sample | Measure | Add sample | Incubate | Measure |
|--|---|---|---|---|
|  |  |  |  |  |
| Insert a swab used for collecting sample into the assay diluent tube. Mix the swab and wait for 1 minute for the mixture to settle down. | [STANDARD TEST] | [READ ONLY] | Leave the test device for 10 minutes. Note that the device should not be left for more than 13 minutes. | Insert the test device into the V200 Analyzer. |

| | | |
|-----------------|---------------------------|----------------------------|
| Ab Titer | Poor immune status | Protective immunity |
| | < Titer 3 | ≥ Titer 3 |